


Accu-Sheath



SHEATH																																
Model	French Size	Sheath Length	Dilator	Guidewire			H Catheter	Seldinger Needle	Wing	Accu Advance	Thumb Scope																					
				Length (cm)	Type	Diameter																										
Femoral	SWSC 4Fr	8	(Pneumocath type)	45	H (Hawley)	18 (0.018")	J	0	0	0	0	0																				
	SWSC 4.5Fr	10	(Radial type)	70	P (Hawley)	250 (0.025")	S (Straight)	1	1	1	1	1																				
Radial	SWSC 5Fr	12	(Radial type)		S (Straight)	350 (0.035")	C (Curved)																									
	SWSC 5.5Fr	14				38 (0.038")																										
	SWSC 6Fr	16																														
	SWSC 6.5Fr	17																														
	SWSC 7Fr	19																														
	SWSC 7.5Fr	21																														
SWSC 8Fr	23																															
	SWSC 8Fr	25																														
<table><tr><th>Guidewire Diameter</th><th>H Catheter</th><th>Seldinger Needle</th><th>Wing</th></tr><tr><td>18 (0.018")</td><td>205</td><td>210</td><td>Transparent</td></tr><tr><td>250 (0.025")</td><td>205</td><td>205</td><td>Transparent</td></tr><tr><td>350 (0.035")</td><td>195</td><td>185</td><td>Pink</td></tr><tr><td>38 (0.038")</td><td>195</td><td>185</td><td>Pink</td></tr></table>													Guidewire Diameter	H Catheter	Seldinger Needle	Wing	18 (0.018")	205	210	Transparent	250 (0.025")	205	205	Transparent	350 (0.035")	195	185	Pink	38 (0.038")	195	185	Pink
Guidewire Diameter	H Catheter	Seldinger Needle	Wing																													
18 (0.018")	205	210	Transparent																													
250 (0.025")	205	205	Transparent																													
350 (0.035")	195	185	Pink																													
38 (0.038")	195	185	Pink																													
Product Item		Order No. Details (Components)																														
Radial (SWSC4Fr)	08L-4SP18C-01110-14Fr 8cm, Radial Long Tip type dilator, 14cm plastic guidewire 0.018" curved type, w/o lv catheter, w/ Accu-Advance, w/ Wing, w/ Accu-Advance, w/ mass handle																															
Femoral (SWSC8Fr)	25F-70535 18000-8Fr 25cm, Femoral Sumpson type dilator, 70cm spring guidewire 0.035" J type, w/ lv catheter, w/ Accu-Advance, w/ Wing, w/ Accu-Advance, w/ mass handle																															

-Sheath System

