

CTO Balloon

Simpass CTO


# mpass

PTCA BALLOON CATHETER



simpass for **CTO** interventions

shaft  
design &  
material



0,016" Entry Profile  
With Atraumatic Soft Tip

that  
folding


**Designed Tapered Down and Hybrid Polymer Distal Shaft**

In Distal Part, 2,5 Fr Profile For Deliverability

In Proximal Part, 2,7 Fr Pusher Body Segment For Superior Pushability

In Distal Part, Soft Polymer Material For Superlative Flexibility





Conical Shaped Shaft Tip For CTO Lesions

0.016" Entry Profile

Comfortable Lesion Crossing With Single Marker

2.5 Fr Flexible Distal Shaft

Pusher Body Design at Conversion Area of Distal Shaft for Enhanced Pushability

Proximal Shaft That Has Increased Kink Resistance

Wide Range Of Product Sizes

